Erkocryl klar/rosa

35 °C
95 °F
41 °F
5 °C

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EN: Instructions for use

Intended use: Thermoforming sheets for the manufacture of apparatus for dental and orthodontic treatment.

Field of application: Erkocryl is thermoformed to fabricate intra-oral appliances such as orthodontic plates, interim dentures, denture bases, dressing or compression plates. Fabrication, see thermoforming technique brochure and videos, www.erkodent.com  

Contraindications: In case of already known allergic reaction to plastics a specific allergy to Erkocryl (PMMA) should be excluded.

Precautionary measures: Erkocryl is only recommended for the above mentioned devices. Observe storage instructions. Attention, sensitive to humidity: Close the aluminium bag immediately.

Attention: The use is limited to the production of dental and orthodontic appliances. Use is the responsibility of the dentist. Use only on the recommendation or prescription of an authorised dentist. For single use only (thermoforming). Allergic reactions are unlikely, but cannot be excluded. Incorrect processing can lead to cracks, fractures and sharp edges. There is a risk of aspiration in the case of small breakages.

Erkocryl (approx. values): density 1.17 g/cm³, modulus of elasticity 2400 MPa, elongation at break 22%, water absorption 0.3%, temperature resistance 96 °C

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