Erkolen

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EN: Instructions for use

Intended use: Thermoforming sheets for the manufacture of apparatus for dental and orthodontic treatment.

Field of application: Erkolen is thermoformed to fabricate casting objects, as moulded pieces for temporary appliances, provisional crowns and bridges, bracket transfer/etching masks as well as spacer for fluoride gel.

Fabrication, see thermoforming technique brochure and videos, www.erkodent.com  

Contraindications: In case of already known allergic reaction to plastics a specific allergy to Erkolen (PE) should be excluded.

Precautionary measures: Erkolen is only recommended for the above mentioned devices. Observe storage instructions.

Attention: The use is limited to the production of dental and orthodontic appliances. Use is the responsibility of the dentist. Use only on the recommendation or prescription of an authorised dentist. For single use only (thermoforming). Allergic reactions are unlikely, but cannot be excluded. Incorrect processing can lead to cracks, fractures and sharp edges. There is a risk of aspiration in the case of small breakages.

Erkolen (approx. values): density 0.92 g/cm³, modulus of elasticity 260 MPa, elongation at break 450 %, water absorption 0.1 %, temperature resistance 68 °C

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