# Kneton-Lab 95

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
# Kneton-Lab 95

Laboratory moulding material on polysiloxan basis, condensation cured, very high consistency, mouldable

Kneadable silicone putty to form sectional moulds for dental repairs in the laboratory, suitable when using cold and hot curing acrylic resins, especially suitable for the fabrication of model bases in case of repairs of breakage, etc.


Consult accompanying documents

Keep away from sunlight


Read instructions for use

Recommended storage temperature

Keep dry

Recommended working temperature

## Dosage:

To 1 level measuring spoon of **Kneton-Lab 95** (17 g) 2.0 cm of hardener paste (0.5 g) are added.

## Application:

Mix **Kneton-Lab 95** with the hardener paste for 30 - 45 sec.

## Tip:

Higher room temperatures and more hardener will accelerate the setting.

## Attention:

Avoid eye contact with the hardener;  
Avoid stains in clothes.

## Technical data:

- mixing can be effected manually or by machine
- mixing time **30 – 45 sec.**
- total working time **1.5 min.**
- setting time **3 min.**
- hardness Shore-A **app. 95**

The user of our products is the only person responsible for the use and the application of our products. Our verbal and written advice has to be considered as advice without obligation also regarding possible protection rights of third parties.

The user is not from the obligation to check our products personally concerning suitability for the intended application. The perfect quality of our products is guaranteed according to our General Terms and Conditions.

In any case and for all damage our liability to the value of the goods supplied by us.