

## Product datasheet Blocking out wax



### 1. Manufacturer information

Trade name:	Blocking out wax
Intended use:	Wax for blocking out
Manufacturer:	Erkodent Erich Kopp GmbH Siemensstraße 3 72285 Pfalzgrafenweiler Germany Tel.: +49 7445 8501-0

### 2. Intended use

Special wax with very high melting point for the thermoforming technique to block out plaster bubbles.

### 3. Composition

CAS-No.:	63231-60-7
Designation:	microcrystalline hard wax

### 4. Properties

Technical properties:

Properties	Guideline	Value
Form	ISO 4823:2015	solid
Colour	-	lilac: titanoxide transparent: oxid free
Odour	-	characteristic
Density	ISO 1183	0,93 g/cm <sup>3</sup>
Melting point	-	70 – 110°C
Viscosity	-	< 30 mPas
Minimum durability	-	36 months

### 5. General information

Storage instructions:

Keep away from sunlight. Keep dry.  
Recommended storage temperature: 5°C – 35°C

The information on this product data sheet corresponds to the best of our knowledge at the time of printing. The information is not guaranteed assurance of product properties and does not constitute a contractual legal understanding. Input errors and mistakes reserved.