

Product datasheet Erkogum



1. Manufacturer information

Trade name:	Erkogum
Intended use:	Blocking out material
Manufacturer:	Erkodent Erich Kopp GmbH Siemensstraße 3 72285 Pfalzgrafenweiler Germany Tel.: +49 7445 8501-0

2. Intended use

Permanently kneadable, reusable blocking out material that withstands the heat during thermoforming. For the protection of residual teeth and to block out undercuts.

3. Properties

Technical properties:

Properties	Guideline	Value
Form	ISO 4823:2015	paste viscous
Colour	-	transparent, lilac
Odour	-	inodorous
Density	ISO 1183	1,2 g/cm ³
Decomposition temperature	-	> 250°C

4. General information

Storage instructions:

Keep away from sunlight. Keep dry.

Recommended storage temperature: 5°C – 35°C

Recommended working temperature: 15°C – 25°C

The information on this product data sheet corresponds to the best of our knowledge at the time of printing. The information is not guaranteed assurance of product properties and does not constitute a contractual legal understanding. Input errors and mistakes reserved.