

Product datasheet Erkoskin



1. Manufacturer information

Trade name:	Erkoskin
Intended use:	Protective film
Manufacturer:	Erkodent Erich Kopp GmbH Siemensstraße 3 72285 Pfalzgrafenweiler Germany Tel.: +49 7445 8501-0

2. Intended use

As spacer when fabricating bleaching trays. To relieve tension of the gingival margin in the thermoforming technique. As protection of plaster dies against trim water, high-glass metal parts when sandblasting and high-gloss ceramic parts when etching.

3. Properties

Technical properties:

Properties	Guideline	Value
Form	ISO 4823:2015	liquid
Colour	-	lilac
Odour	-	Ammonia (NH ₃)
Density	ISO 1183	0,947 g/cm ³
pH	-	10 – 12
Minimum durability	-	36 months

4. General information

Storage instructions:

Keep away from sunlight. Keep dry. Put lids on containers immediately after use.

Recommended storage temperature: 5°C – 35°C

Recommended working temperature: 15°C – 25°C

The information on this product data sheet corresponds to the best of our knowledge at the time of printing. The information is not guaranteed assurance of product properties and does not constitute a contractual legal understanding. Input errors and mistakes reserved.