

## Product datasheet Sitran-F



### 1. Manufacturer information

Trade name:	Sitran-F
Intended use:	Correction material for the two-stage impression technique
Manufacturer:	Erkodent Erich Kopp GmbH Siemensstraße 3 72285 Pfalzgrafenweiler Germany Tel.: +49 7445 8501-0

### 2. Intended use

Condensation-cured high precise impression material on polysiloxan basis for the fabrication of high precise correction and two-stage impressions.

### 3. Properties

#### Technical properties:

Properties	Guideline	Value
Form	ISO 4823:2015	paste viscous
Colour	-	blue
Odour	-	characteristic
Density	ISO 1183	1,23 g/cm <sup>3</sup>
Minimum durability	-	24 months

### 4. General information

#### Storage instructions:

Keep away from sunlight. Keep dry.  
Recommended storage temperature: 5°C – 35°C  
Recommended working temperature: 15°C – 25°C

#### Instructions for cleaning and maintenance:

Do not bring into contact with disinfectants or similar.

The information on this product data sheet corresponds to the best of our knowledge at the time of printing. The information is not guaranteed assurance of product properties and does not constitute a contractual legal understanding. Input errors and mistakes reserved.