

All the performance you expect

When comparing a Dentist fabricated Playsafe mouthguard with shop brought boil & bite mouthguards the differences are easy to see and easier to feel. A Playsafe mouthguard fits exactly onto your teeth, doesn't move and allows clear breathing and speech.

Boil & bite mouthguards can never fit as accurately or be as comfortable as a professionally fitted Playsafe mouthguard and give a false sense of security by lacking the ability to dissipate the force of serious sporting impacts, the time when it counts the most.





Playsafe Mouthguard



Boil & bite mouthguard

Independently tested on and off the field, a Playsafe mouthguard uses proven and innovative design to dramatically reduce your risk of sports dental injuries, giving you the confidence to compete at your best




Your Sport, your Safety, your Playsafe



Playsafe has a mouthguard suited for every player from juniors to sports competitors at the highest levels. The range of Playsafe mouthguards lets you select the mouthguard that is exactly right for you. You never need to purchase a mouthguard that is too bulky, uncomfortable or does not give you the protection you expect. Generous Private Health Fund rebates generally apply for Playsafe mouthguards, ask your Dentist

Available from



© Erkodent Australia 2018



PLAYSAFE
COMPETITION MOUTHGUARDS



THE ULTIMATE SPORTS PROTECTION

PLAYSAFE™

COMPETITION MOUTHGUARDS

The perfect mouthguard for every sport

If you would like to order using this form please:

1. Select your mouthguard type from below
2. Choose your colour(s) and/or graphics
3. Complete your details
4. Return the completed form to your Dentist

For high impact and professional sports

- ☐ HEAVY PRO ☐ TRIPLE
☐ LIGHT PRO ☐ TRIPLE LIGHT


Made from 3 layers, using a hard internal layer. For sports where extremely hard or pointed impacts may be received. The hard layer in the Playsafe Heavy and Light Pro extends over the outer teeth area only while the Triple and Triple Light incorporates the hard layer over the entire mouthguard, this gives a slightly firmer feel but means that there can be far less bulk in the tongue area (see diagram) These mouthguards are best suited to stable adult teeth, use the Light versions if space or bulk is a concern. For sports such as Senior Rugby, Hockey, AFL, Ice Hockey, Karate, Squash, Kickboxing etc

For club sports and junior players

- ☐ MEDIUM ☐ LIGHT

Playsafe Medium and Light are dual laminated mouthguards. The Medium is designed for general use and sports with likely impacts from objects with a larger surface area eg Boxing, Basketball, Soccer, Judo, Wrestling etc. Choose the Playsafe Light for children (up to ~13yrs) and sports that require a mouthguard to protect lower teeth from damage (eg motor-sports). Flexible and comfortable and can be readapted to changing teeth.

Colours




Freestyle



Multi-colours

Create your mouthguard with 2, 3 or 4 colours. (N/A for Playsafe Triple) Indicate your choice below from colours 1-16



Playsafe 4u Graphics

Add graphics to your mouthguard. Mark the position and code # on the mouthguard outline below. Our graphics are updated regularly so ask your Dentist to see the current selection.



#16



#03



#04



#05



#18



#01 white



#02 black



#17



#06



#07



#08



#10



#11



#12



#13



#14



#09



#15



Your details

If you wish to provide your contact details they will be permanently sealed in your Playsafe mouthguard at no additional cost

Name

Telephone